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

Mask component	Material
Headgear	Elastane, polyester, PU (polyurethane), UBL loop, cotton, PA (polyamide)
Anti-asphyxia valve	TPE, PP (polypropylene)
Ripcord	POM, polyester, PA (polyamide)

All mask components are free from latex, PVC (polyvinyl chloride) and DEHP (diethylhexyl phthalate).

11 Warranty

Weinmann gives the customer a limited manufacturer warranty on new genuine Weinmann products and any replacement part fitted by Weinmann in accordance with the warranty conditions applicable to the product in question and in accordance with the warranty periods from date of purchase as listed below. The warranty conditions can be downloaded from www.weinmann-medical.com on the Internet. We can also send you the warranty conditions on request. Please note that any claim under warranty and liability will be void if neither the accessories nor original replacement parts recommended in the instructions for use are used.

In the event of a claim under warranty, please contact your specialist dealer.

Product	Period of guarantee
Masks including accessories	6 months

12 Declaration of conformity

Weinmann Geräte für Medizin GmbH + Co. KG (Kronsaaßweg 40, 22525 Hamburg, Germany) hereby declares that the product complies with the relevant regulations of Directive 93/42/EEC governing medical devices. The complete text of the declaration of conformity can be found at: www.weinmann-medical.com

1 Introduction

1.1 Intended use

The JOYCEone Full Face masks are used for treating sleep apnea and for non-invasive and non-life-sustaining ventilation of patients with respiratory insufficiency. They serve as a connecting element between the patient and the therapy device. Only use non-vented full face masks in combination with therapy devices equipped with an active exhalation valve.

1.2 Contraindications

In the case that any of the following symptoms are present the mask may not be used, or only with particular caution:

erosions and ulcerations, skin allergies, flush of the facial skin, pressure marks on the face, claustrophobia, anxiety, deformities of the face or nasopharynx, intake of medication which may induce vomiting, the necessity for immediate intubation. Please also observe the contraindications in the instructions for use of your therapy device.

1.3 Side effects

Nasal congestion, dry nose, morning xerostomia (dry mouth), feeling of pressure in the sinuses, irritated conjunctiva, skin rashes, pressure marks on the face, irritating noises when breathing.

2 Safety

Please read these instructions carefully. They are a constituent part of the mask and must be available at all times.

Risk of injury from the supply of oxygen!

Oxygen can deposit in clothing, bed linen, and hair. In conjunction with smoking, naked flames, or electrical equipment, it can cause fires and explosions.

⇒ Do not smoke.

⇒ Keep away from naked flames.

Risk of injury from re-inhaling CO₂!

If the mask is used incorrectly, CO₂ may be re-inhaled. ⇒ For vented versions: Do not seal off the exhalation system.

⇒ Only apply the mask for a prolonged period if the therapy device is running.

⇒ Only use the mask within the quoted therapy pressure range.

⇒ Do not use the mask on patients who cannot take it off themselves.

⇒ Only use non-vented full face masks in combination with therapy devices equipped with an active exhalation valve.

Risk of injury from slipping of the mask!

If the mask slips or falls off, the therapy is ineffective. ⇒ Monitor patients with restricted spontaneous respiration.

⇒ Activate low pressure/leakage alarms on the therapy device.

Risk of injury from anesthetic gases!

Anesthetic gases can escape through the exhalation valve and put third parties at risk.

⇒ Do not use the mask during anesthesia.

3 Product description

3.1 Overview

Illustrations of the individual components can be found on the front cover.

Key

1. Headgear
2. Forehead cushion
3. Forehead support
4. Retaining ring
5. Elbow
6. Anti-asphyxia valve (vented only)
7. Plug for pressure measurement port (vented only)
8. Pressure measurement port (vented only)
9. Rotating sleeve (vented only)
10. Mask body
11. Headgear clip
12. Mask cushion

3.2 Compatible devices

With some combinations of devices, the actual pressure in the mask will not correspond to the prescribed therapy pressure, even if the therapy device is displaying the correct pressure. Have the combination of devices adjusted by a physician or specialist dealer so that the actual pressure in the mask corresponds to therapy pressure.

3.3 Exhalation system

Vented versions

The vented versions (transparent elbow and rotating sleeve) have an integrated exhalation system. The retaining ring and mask body are shaped in such a way that a gap forms between these parts. The exhaled air can escape through this gap.

Non vented versions

The non vented versions (blue retaining ring and elbow) do not have an exhalation system. Non vented full face masks may only be used in combination with therapy devices equipped with an active exhalation valve as well as alarms and safety systems in case the device should fail. If you use the non vented version with an external exhalation valve, please observe the corresponding instructions for use.

3.4 Anti-asphyxia valve

▲ WARNING

Risk of asphyxia if anti-asphyxia valve not working properly!

Residues can block the valve, resulting in re-inhalation of CO₂.
⇒ Check that the openings of the anti-asphyxia valve are clear before every use.

Should the therapy device fail, the anti-asphyxia valve opens so that the patient can breathe ambient air.

3.5 Pressure measurement port

The mask has a pressure measurement port which can be used for measuring therapy pressure or supplying oxygen. If you are not using this, close it with the plug for pressure measurement port.

4 Operation

The following figures show you how to apply, adjust, remove, dismantle and assemble the mask:

Figure 1: Putting on the mask

Figure 2: Adjusting the mask

Figure 3: Removing the mask

Figure 4: Dismantling the mask

Figure 5: Assembling the mask

5 Hygienic preparation

▲ WARNING

Risk of injury from inadequate cleaning!

Residues may congest the mask, impair the integrated exhalation system, and jeopardize therapy success.

⇒ For patients with a compromised immune system or particular background of illness, disinfect mask components daily following consultation with the physician.

5.1 Cleaning the mask

- Dismantle the mask (see figure 4).
- Clean mask in accordance with the table below:

Action	Daily	Weekly
Wash mask components with warm water and a mild detergent.	X	
When washing mask parts (exception: mask cushion and forehead cushion), clean thoroughly with a cloth or a soft brush.		X
Wash headgear by hand.		X

- Rinse all parts with clear water.

- Allow all parts to air-dry.

- Perform a visual inspection.

- If necessary: Replace damaged parts.

- Assemble the mask (see figure 5).

i

Discolorations of mask components do not impair the functionality of the mask.

5.2 Change of patient

You must subject the mask to a hygienic preparation in the event of a change of patient. For information on hygienic preparation in the event of a change of patient, please refer to our Internet brochure at www.weinmann-medical.com. If you do not have Internet access, please contact Weinmann.

6 Disposal

You can dispose of all parts in domestic waste.

7 Faults

Fault	Cause	Remedy
Facial pressure pain.	Mask fits too tightly.	Loosen headgear.
Draft in the eyes.	Mask fits too loosely. Mask does not fit.	Tighten headgear. Contact your authorized dealer.

Quoted two-figure noise emission value according to ISO 4871 (vented only):	15 dB(A) 23 dB(A) 3 dB(A)
- Sound pressure level	
- Uncertainty factor	
Service life	Up to 12 months
Standards applied	EN ISO 17510-2:2009


The materials of the mask deteriorate if exposed to e.g., aggressive detergents. In individual cases, therefore, it may be necessary to replace mask parts sooner.

CE 0197

The right to make design modifications is reserved.

9 Pressure/flow characteristic curve

The pressure/flow characteristic curve shows the leak flow as a function of therapy pressure.

10 Materials

Mask component	Material
Mask cushion, forehead cushion	Silicone, silicone coating
Plug for pressure measurement port, spring of forehead support	Silicone
Retaining ring (hard component)	PP (polypropylene)
Retaining ring (soft component)	TPE (thermoplastic elastomer)
Headgear clip	PA (polyamide), POM (polyoxymethylene)
Forehead support, rotating sleeve, mask body, elbow	PA (polyamide), PC (polycarbonate)

Fault	Cause	Remedy
	Mask is not adjusted correctly.	Re-adjust mask (see figure 2).
	Mask cushion is damaged.	Replace mask cushion.
	Tube system is leaking.	Check connection and fit of tubes.
Therapy pressure is not reached.	Air is escaping from pressure measurement port (vented only).	Close pressure measurement port with plug (included in scope of supply).
	Anti-asphyxia valve incorrectly fitted (vented only).	Fit anti-asphyxia valve correctly (see figure 5).
	Anti-asphyxia valve defective (vented only).	Replace anti-asphyxia valve.

8 Technical data

Product class according to 93/42/EEC	IIa
Dimensions (W x H x D)	120 mm x 150 mm x 100 mm
Weight Vented	100 g
Non vented	95 g
Therapy pressure range	4 hPa - 30 hPa
Hose connection: Cone according to EN ISO 5356-1 Vented	Ø 22 mm (male) Ø 22 mm (female)
Non vented	
Pressure measurement port	Ø 4 mm
Width of headgear clips	24 mm
Temperature range:	
Operation	+5°C to + 40°C
Storage	-20°C to + 70°C
Flow resistance vented at 50 l/min	0.05 hPa
at 100 l/min	0.28 hPa
Flow resistance non vented at 50 l/min	0.08 hPa
at 100 l/min	0.32 hPa
Anti-asphyxia valve flow resistance (vented only)	
Inspiration at 50 l/min:	0.5 hPa
Exhalation at 50 l/min:	0.5 hPa
Switching pressure for anti-asphyxia valve (vented only)	
• Open:	1 hPa
• Close:	2 hPa

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